

SP 200 with 12-Slot Collating Control Center

Robotic Prescription Dispensing System

ScriptPro's SP 200 with 12-Slot Collating Control Center is a proven and reliable prescription dispensing system that automates the repetitive, manual dispensing tasks most subject to human error and automatically collates all prescriptions for each patient.

SP 200 interfaces with the pharmacy computer system to fill, label, and collate up to 150 prescriptions per hour. Each slot displays the patient name for easy vial retrieval.

The system can be configured to collate the prescriptions in various ways, such as for a given patient or dispensing priority. Electronic displays on both sides of the control center show the contents of the slot.

SP 200 contains 200 universal dispensing cells, which are easily calibrated on-site by pharmacy staff, three vial dispensers or two bulk load dispensers, and the 12-Slot Collating Control Center. It dispenses tablets and capsules of all shapes and sizes into standard pharmacy vials. There is no drug cross-contamination because the system fills directly from the dispensing cell into the vial. The system also prints and applies the prescription and auxiliary labels and collates the uncapped vials into slots for final inspection using on-screen drug image verification.

Barcodes are used throughout the process for accuracy and quality control. Pharmacy staff manages the dispensing process from prescription entry to approval of the finished product.



SP 200 with 12-Slot Collating Control Center and Bulk Load Vial Dispensers

**3 Vial Dispensers
(optional)**


ScriptPro®

Pharmacy Automation


For more information call: 800.606.7628 • 913.384.1008 • www.scriptpro.com

SP 200 with 12-Slot Collating Control Center


Top View



Front View



Back View



Specifications*

Weight:

- A. Robotics Cabinet: 850.5 lbs (385.8 kg)
 - 1. Vial Dispensers: 166.5 lbs (75.5 kg)
 - 2. Bulk Load Dispensers: 87.5 lbs (39.7 kg)
 - 3. Spiral Vial Dispenser: 64 lbs (29 kg)
- B. Electrical Cabinet: 326 lbs (147.9 kg)
- C. Collating Control Center: 289 lbs (131 kg)

Robot Footprint: 24.08 sq. ft. (2.24m²)

Electrical Input: 100-127 VAC, 50/60 Hz, 12 AMP (robot only; includes UPS), 16 AMP with all peripherals, 1 Phase, Uninterruptible Power Supply included. Dedicated 20 AMP power outlet located within 10 feet (3m) of power connection.

Communications: Serial, parallel, or ethernet pharmacy computer connection. Firewall provided for network connectivity.

Remote Support Requirements:

1. VPN over high-speed Internet connection. Also requires access to backup telephone line as specified in #2 below, such as a fax line.

OR

2. Modem-Dedicated analog RJ-11 telephone outlet accessible to outside line and capable of sending and receiving calls.

Barcode Scanner - Caution: Class IIa Laser Product. Do not stare into beam. 650nm visible laser light diode - .663 mW maximum output.

Computer/Software: Client/Server system on ScriptPro-installed network; Intel Core 2 Duo PCs with dual SATA drives. TCP/IP socket connections; USB connections; Software built with Visual C++; Windows XP/MSDE. Redundant data and transaction log backups. Customer access limited to SP Central user interface.

Monitor: 17" touchscreen with articulating arm

Temperatures:

- Operating: 41°F to 86°F (5°C to 30°C)
- Storage: 5°F to 104°F (-15°C to 40°C)
- Relative Humidity: 20-85%
- Altitude: Up to 9,843 feet (3000m)
- Dissipated Heat: 1,926 BTU/HR max

Cell Capacity:

- Single Cell: 676 cc/183 dram
- Extended Cell: 1900 cc/514 dram

Peripheral Products Supported:

See product sheet on peripheral products

System has 12 slots with a maximum capacity of 7 vials per slot.

ANSI Barcode Print Quality Grade: B

Sound Level: 77 db at front edge of Control Center

Training: Up to four Full Time Equivalent (FTE) training days included with original lease/purchase. Approx. two FTE days for Training and Setup and approx. two FTE days for Go-Live and Support.

One "FTE" day is defined as one trainer on-site for one day for up to 10 hours. If more than one trainer is on-site, it is counted as multiple FTE days.

* Specifications subject to change without notice.

** Requires one additional inch clearance for installation.



ScriptPro

Pharmacy Automation

For more information call: 800.606.7628 • 913.384.1008 • www.scriptpro.com

November 2011